

Patiëntenlijst

On this page you can view a list of patients and have access to their pictures and reports!

With this button you can call **help** to this site!

WEB access > Beeld- en verslaginzage > Patiëntenlijst

Enter time range!

Click on this button to open the calendar function. Using the calendar you can easily select a date, which is then automatically inserted in the corresponding text field.

Select filter parameters

Search value to the filter

Vanaf: tot: Filter: Waarde:

As soon as you have filled in all the desired search criteria, click on this button to **load the list!**

Click on one of the links on the patient (not report and image) to get to the patient page.

This contains all the important information about the patient. You can also view images and manage and play dictations.

As soon as a report or image is available, there are links to each patient. By pressing the respective link, you can retrieve the data for the patient.

in de patiënt, om de patiënt informatie weer te geven
in XR XRWEB acces of XR PACS imagepro (ActiveX) openen
als u een windows 8 bezitten, bitte das ImagePro-Setup (vom 12.01.2015) neu runterladen, mit Hilfe des neuen Setup die alte Version deinstallieren und die neue Version installieren. Folgen Sie der Anleitung Hilfe zu ActiveX bzw. Algemene hulp naar beeldinzage
Zip-File: beelden als Zip downloaden
OsiriX:URLs-File downloaden, die in OsiriX geopend kan worden

Vorname	Geboren	Verz	Verwijzer	Onderzoek	Onderzoeking	Verslag	Beeld
L	11.03.1957		[redacted]	07.04.2015	cervicale wk	verslag	geen
R	12.06.1993	Atrium	[redacted]	07.04.2015	POLS/HAND/MIN R	uitgever	gearchiveerd Z OsiriX
GAM	15.03.1953		[redacted]	07.04.2015	schouder/humerus r, echo/schouder r.	verslag	geen
B	20.05.1968	Atrium	[redacted]	07.04.2015	ECHO/SCHOUDER R., SCHOUDER/HUMERUS R, CERVICALE WK	klaar	gearchiveerd Z OsiriX
Heer	[redacted]	J	[redacted]	07.04.2015	thorax	verslag	geen



durchschaubare IT für die Radiologie

